



Fentanyl Confirmation in Oral Fluids by Laser Diode Thermal Desorption (LDTD) – MS/MS

Serge Auger¹, Alex Birsan¹ & Dean Fritch²

¹Phytronix Technologies, Quebec, Canada, ²OraSure Technologies, Bethlehem, USA

Keywords: Fentanyl Analysis, Oral Fluids, Ultra High-Throughput, LDTD

Introduction

Drug testing in Oral Fluids is a constantly evolving analysis procedure which benefits from increasingly sensitive methods of detection. Testing for drugs of abuse in oral fluids can strongly benefit the criminal justice field as a less invasive and cost-effective approach for drug detection when compared to blood or urine sampling. In a clinical environment, oral fluids can be used in patient screening for rapid confirmation of the presence or absence of orally administered drugs.

The LDTD ion source uses an infrared laser diode to desorb samples that have been dried onto a 96-well LazWell™ plate. The rapid desorption produces neutral species which are carried into a corona discharge region to undergo an efficient protonation and are subsequently transferred directly into the mass spectrometer for detection.

Oral Fluid Collection

- The Intercept® device by OraSure is used for saliva collection. Standard curves and QC's are prepared in the Oral Fluid Calibration Buffer.



Figure 1: Intercept® Oral Fluid Drug Test

LDTD-MS/MS System



Figure 2: LDTD system on Thermo Vantage Mass Spectrometer.

Sample Method

Extraction Procedure

- 100 µL Oral Fluid Calibration Buffer
- 20 µL IS (Fentanyl-d5 at 20 ng/mL in MeOH)
- 100 µL NaOH (0.1N in Water)
 - Mix
- 600 µL Ethyl Acetate*
 - Mix and centrifuge (2 min. / 14000 rpm)
- Spot 5 µL of organic phase in LazWell plate
 - Evaporate to dryness

*Organic phase can be evaporated and reconstituted to further concentrate the sample

LDTD-MS/MS Parameters

LDTD

Gas Flow:	3 L/min	
Laser pattern:	Time (s)	Power (%)
	0	0
	2	0
	5	45
	7	45
	7.1	0
	8	0

MS/MS Method

	Transition	CE	S-Lens
Fentanyl	337->188	22	120
Fentanyl-d5	342->188	22	120
Mode:	Positive		

Results and Discussion

Linearity Results

As shown in **Figure 3**, excellent linearity ($r^2 > 0.99$) with no signs of carryover effect is achieved in the quantification range (0.5 to 50 ng/mL).




Figure 3: Fentanyl standard curve

	r^2	Slope (ratio area / concentration)	y- Intercept
Run 1	0.9992	0.2946	0.0171
Run 2	0.9969	0.2947	0.0186
Run 3	0.9990	0.2994	0.0307

Table 1: Calibration Curve Parameters

Accuracy and Precision

As shown on **Table 2 and 3**, the inter-run and intra-run accuracy and precision are between 90.1 to 106.8% and 0.6 to 7.7%, respectively.

	QC-Low	QC-Med	QC-High
Conc. (ng/mL)	1.25	5	25
N	9	9	9
Mean (ng/mL)	1.29	5.10	25.03
%RSD	7.7	5.8	2.3
%Nom	103.3	102.0	100.1

Table 2: Inter-run precision and accuracy for Fentanyl

	LLOQ	QC-Low	QC-Med	QC-High	ULOQ
Conc. (ng/mL)	0.5	1.25	5	25	50
N	3	3	3	3	3
Mean (ng/mL)	0.45	1.34	5.02	25.21	49.17
%RSD	4.2	2.3	3.9	0.9	0.6
%Nom	90.1	106.8	100.4	100.8	98.3

Table 3: Intra-run precision and accuracy for Fentanyl

Detection limit (LOD)

A detection limit of 0.25 ng/mL can be reached with a blank interference of 20.0% at this concentration.

Stability Verification

Following the extraction process, all samples were stored at 4°C to evaluate the wet stability of the drug. After 42h40, all samples were re-spotted and analyzed. Linearity, precision and accuracy were evaluated to determine the stability. **Table 4** shows that a wet stability of 42h40 is obtained with good precision and accuracy of LOQ standard.

The stability of dry samples in LazWell plate was also determined. All standards and QCs are spotted, dried and kept at room temperature for 42h20. Then, standards and QCs were analyzed and the linearity, precision and accuracy are verified. **Table 4** shows the dry stability results and the storage conditions of the LazWell.

	Wet Stability	Dry in LazWell (RT)
Time (h)	42h40	42h20
Temp. (°C)	4°C	RT
Conc. (ng/mL)	1	1
N	3	3
Mean (ng/mL)	0.46	0.50
%RSD	17.9	11.0
%Nom	92.4	99.5

Table 4: Stability Results for Fentanyl

Correction Factor

Values reported represent diluted oral fluid. To convert to whole saliva, you must multiply by a factor of 3X.

Conclusions

The ease of use of the Intercept® oral fluid sampling device from OraSure provides an accurate and fast sampling method for many drugs of abuse. The combination of the oral fluid extraction procedure with the analysis speed of the LDTD-MS/MS is an ideal solution in high-throughput drug analysis.

A fast, sensitive and reproducible method for the analysis of fentanyl in oral fluid matrix is achieved using a simple buffer extraction method combined to the speed of analysis of the LDTD-MS/MS with a total sample-to-sample analysis time of **8 seconds**.