

# Solid-Phase Extraction of MDMA, MDA and MDEA in Urine and analysis by LDTD-MS/MS

Serge Auger<sup>1</sup>, Alex Birsan<sup>1</sup>, Annick Dion-Fortier<sup>1</sup> & David Dubé<sup>2</sup>

<sup>1</sup>Phytronix Technologies, Québec, Canada, <sup>2</sup>SiliCycle, Québec, Canada

**Keywords:** MDMA, MDA, MDEA, Urine, Solid-Phase Extraction, LDTD

## Introduction

Analysis of certain drugs of abuse in urine can require a sample clean-up step to reduce the interference effect from the matrix. To obtain an optimal sample clean-up, the SiliaPrep™ CleanDrug SPE cartridges are used in the extraction procedure prior to the Laser Diode Thermal Desorption (LDTD) analysis.

The LDTD ion source uses an infrared laser diode to desorb samples that have been previously dried onto a 96-well LazWell™ plate after sample preparation extraction. The rapid desorption produces neutral species which are carried into a corona discharge region to undergo an efficient protonation and are subsequently transferred directly into the mass spectrometer for detection.

## Solid Phase Cartridge

The SiliaPrep CleanDrug cartridge is used for the sample extraction procedure.



Figure 1: SiliaPrep CleanDrug SPE cartridge

SiliaPrep CleanDRUG Formats		
Formats	Qty / Pk	Product number
1 mL / 50 mg	100	SPEC-R651230B-01B
1 mL / 100 mg	100	SPEC-R651230B-01C
3 mL / 200 mg	50	SPEC-R651230B-03G
3 mL / 500 mg	50	SPEC-R651230B-03P
6 mL / 500 mg	50	SPEC-R651230B-06P
6 mL / 1g	50	SPEC-R651230B-06S
2 mL / 50 mg	1	96W-R651230B-B
2 mL / 100 mg	1	96W-R651230B-C

Table 1: SiliaPrep CleanDRUG product number

## LDTD-MS/MS System



Figure 2: LDTD system on AB SCIEX 5500 Qtrap Mass Spectrometer

## Sample Method

### Extraction procedure

**Cartridge:** SiliaPrep CleanDrug (1 mL / 100 mg)  
**Activation:** 1 mL MeOH  
 1 mL Water  
 1 mL Na Acetate (100 mM, pH 6) in Water  
**Load:** 200 µL sample  
 40 µL IS (MDMA-d5 at 200 ng/mL in MeOH)  
 600 µL Na Acetate (100 mM, pH 6) in Water  
**Wash 1:** 1 mL Water  
**Wash 2:** 1 mL MeOH  
**Elution:** 1 mL EtAc/IPA/NH<sub>4</sub>OH (78/20/2)  
 After elution, add 40 µL Formic Acid. Mix\*  
 Spot: 2 µL in LazWell plate

\*Organic phase can be evaporated and reconstituted to further concentrate the sample

## LDTD-MS/MS Parameters

### LDTD

Gas Flow:	3 L/min	
Laser pattern:	Time (s)	Power (%)
	0	0
	2	0
	5	45
	7	45
	7.1	0
	8	0

### MS/MS Method

	Transition	CE	DP
MDMA	194->163	12	30
MDMA-d5	199->165	12	30
MDA	180->133	20	30
MDEA	208->163	12	30
Mode:	Positive		

## Results and Discussion

### Linearity Results

As shown in **Figure 3**, excellent linearity ( $r^2 > 0.99$ ) with no signs of carryover effect is achieved within the quantification range (50 to 5,000 ng/ml).




Figure 3: Typical standard curve

### Accuracy and Precision

As shown on **Table 2 and 3**, the inter-run and intra-run accuracy and precision are between 96.3 to 111.5% and 1.8 to 11.8% for all three drugs.

		Conc. (ng/mL)	N	Mean (ng/mL)	%RSD	%Nom
MDMA	QC-Low	125	9	126.0	9.3	100.8
	QC-Med	500	9	491.9	7.1	98.4
	QC-High	2500	9	2502.2	5.1	100.1
MDA	QC-Low	125	9	131.7	6.8	105.3
	QC-Med	500	9	538.4	5.5	107.7
	QC-High	2500	9	2550.0	7.7	102.0
MDEA	QC-Low	125	9	128.6	8.9	102.8
	QC-Med	500	9	531.0	8.0	106.2
	QC-High	2500	9	2565.6	4.9	102.6

Table 2: Inter-run precision and accuracy

		Conc. (ng/mL)	N	Mean (ng/mL)	%RSD	%Nom
MDMA	L L O Q	50	3	49.7	7.8	99.5
	QC-Low	125	3	121.7	8.1	97.3
	QC-Med	500	3	502.7	9.7	100.5
	QC-High	2500	3	2450.0	4.6	98.0
	U L O Q	5000	3	5003.3	2.5	100.1
MDA	L L O Q	50	3	51.0	4.7	102.1
	QC-Low	125	3	124.3	9.5	99.5
	QC-Med	500	3	553.0	1.8	110.6
	QC-High	2500	3	2600.0	10.2	104.0
	U L O Q	5000	3	5013.3	8.3	100.3

		Conc. (ng/mL)	N	Mean (ng/mL)	%RSD	%Nom
MDEA	L L O Q	50	3	48.1	11.8	96.3
	QC-Low	125	3	123.0	7.1	98.4
	QC-Med	500	3	557.7	8.0	111.5
	QC-High	2500	3	2646.7	4.5	105.9
	U L O Q	5000	3	4896.7	4.5	97.9

Table 3: Intra-run precision and accuracy

### Recovery

Recovery at 5,000 ng/mL of concentration for each drug is reported in **Table 4** (N=3).

	MDMA	MDEA	MDA
Recovery (%)	93	100	90

Table 4: Recovery results for all drugs

### Stability Verification

Following the SPE extraction process, all samples were stored at 4°C to evaluate the wet stability of the drugs. After 96h, all samples were re-spotted and analyzed. Linearity, precision and accuracy were evaluated to determine the stability. **Table 5** shows that a wet stability of 96h is obtained with good precision and accuracy of LOQ standard.

The stability of dry samples in LazWell plate was also determined. All standards and QCs are spotted, dried and kept at room temperature for 24h. Then, standards and QCs were analyzed and the linearity, precision and accuracy are verified. **Table 5** shows the dry stability results and the storage conditions of the LazWell.

Time (h)	Wet Stability			Dry in LazWell (RT)		
	MDA	MDEA	MDMA	MDA	MDEA	MDMA
96						
24						
Temp. (°C)	4°C			RT		
Conc. (ng/mL)	50			50		
N	3			3		
Drug	MDA	MDEA	MDMA	MDA	MDEA	MDMA
Mean (ng/mL)	50.6	46.8	49.2	44.7	45.7	53.2
%RSD	16.1	14.2	15.6	1.2	8.6	9.2
%Nom	101.1	93.6	98.4	89.3	91.3	106.5

Table 5: Stability Results for MDEA, MDA and MDMA

### Conclusions

The sample solid extraction using **SiliaPrep CleanDrug** SPE cartridges ensures accurate and precise results with a linear standard curve ( $r^2 > 0.99$ ) for all three drugs.

A fast analysis can be reach using LDTD-MS/MS system. This system allows a total sample-to-sample analysis time of **8 seconds**.