



Butalbital, Secobarbital and Phenobarbital Confirmation in Oral Fluids by Laser Diode Thermal Desorption (LDTD) – MS/MS

Serge Auger¹, Alex Birsan¹ & Dean Fritch²

¹Phytronix Technologies, Quebec, Canada, ²OraSure Technologies, Bethlehem, USA

Keywords: Butalbital, Secobarbital, Phenobarbital, Oral Fluids, Ultra High-Throughput, LDTD

Introduction

Drug testing in Oral Fluids is a constantly evolving analysis procedure which benefits from increasingly sensitive methods of detection. Testing for drugs of abuse in oral fluids can strongly benefit the criminal justice field as a less invasive and cost-effective approach for drug detection when compared to blood or urine sampling. In a clinical environment, oral fluids can be used in patient screening for rapid confirmation of the presence or absence of orally administered drugs.

The LDTD ion source uses an infrared laser diode to desorb samples that have been dried onto a 96-well LazWell™ plate. The rapid desorption produces neutral species which are carried into a corona discharge region to undergo an efficient protonation and are subsequently transferred directly into the mass spectrometer for detection.

Oral Fluid Collection

- The Intercept® device by OraSure is used for saliva collection. Standard curves and QC's are prepared in the Oral Fluid Calibration Buffer.



Figure 1: Intercept® Oral Fluid Drug Test

LDTD-MS/MS System



Figure 2: LDTD system on Xevo-TQMS mass spectrometer

Sample Method

Extraction Procedure

- 100 µL Oral Fluid Calibration Buffer
- 20 µL IS (20 µg/mL in MeOH:Water (75/25))
- 100 µL HCl (0.1N in Water)
 - Mix
- 600 µL Ethyl Acetate*
 - Mix and centrifuge (2 min. / 14000 rpm)
- Spot 5 µL of organic phase in LazWell plate
 - Evaporate to dryness

*Organic phase can be evaporated and reconstituted to further concentrate the sample

LDTD-MS/MS Parameters

LDTD

Gas Flow:	3 L/min	
Laser pattern:	Time (s)	Power (%)
	0	0
	2	0
	5	45
	7	45
	7.1	0
	8	0

MS/MS Method

	Transition	CE	Cone
Butalbital	225->98	30	25
Secobarbital	239->98	20	20
Phenobarbital	233->162	15	25
Mode:	Positive		

Results and Discussion

Linearity Results

As shown in **Figure 3**, excellent linearity ($r^2 > 0.99$) with no signs of carryover effect is achieved within the quantification range (20 to 2,000 ng/mL).




Figure 3: Typical standard curve

Accuracy and Precision

As shown on **Table 1 and 2**, the inter-run and intra-run accuracy is between 88.8 to 112.2% and the precision is between 2.5 to 18.2%.

		Conc. (ng/mL)	N	Mean (ng/mL)	%RSD	%Nom
Butal.	QC-Low	50	9	49.44	5.7	98.9
	QC-Med	500	9	488.10	10.1	97.6
	QC-High	1000	9	991.83	10.1	99.2
Seco.	QC-Low	50	9	48.31	10.4	96.6
	QC-Med	500	9	504.33	8.9	100.9
	QC-High	1000	9	961.94	4.7	96.2
Pheno.	QC-Low	50	9	53.09	7.0	106.2
	QC-Med	500	9	487.14	9.3	97.4
	QC-High	1000	9	937.50	4.5	93.8

Table 1: Inter-run precision and accuracy

		Conc. (ng/mL)	N	Mean (ng/mL)	%RSD	%Nom
Butal.	LLOQ	20	3	19.40	13.2	97.0
	QC-Low	50	3	48.13	5.8	96.3
	QC-Med	500	3	465.50	6.8	93.1
	QC-High	1000	3	960.40	11.6	96.0
	ULOQ	2000	3	2048.00	9.7	102.4
		Conc. (ng/mL)	N	Mean (ng/mL)	%RSD	%Nom
Seco.	LLOQ	20	3	22.43	18.2	112.2
	QC-Low	50	3	48.65	8.3	97.3
	QC-Med	500	3	466.80	9.3	93.4
	QC-High	1000	3	966.73	5.3	96.7
	ULOQ	2000	3	2065.53	6.4	103.3

		Conc. (ng/mL)	N	Mean (ng/mL)	%RSD	%Nom
Pheno.	LLOQ	20	3	21.33	7.6	106.7
	QC-Low	50	3	52.17	6.8	104.3
	QC-Med	500	3	444.13	2.5	88.8
	QC-High	1000	3	939.43	6.4	93.9
	ULOQ	2000	3	2117.37	10.2	105.9

Table 2: Intra-run precision and accuracy

Stability Verification

Following the extraction process, all samples were stored at 4°C to evaluate the wet stability of the drugs. After 73h, all samples were re-spotted and analyzed. Linearity, precision and accuracy were evaluated to determine the stability. **Table 3** shows that a wet stability of 73h is obtained with good precision and accuracy of LOQ standard.

The stability of dry samples in LazWell plate was also determined. All standards and QCs are spotted, dried and kept at room temperature for 71h. Then, standards and QCs were analyzed and the linearity, precision and accuracy are verified. **Table 3** shows the dry stability results and the storage conditions of the LazWell.

	Wet Stability			Dry in LazWell (RT)		
Time (h)	73			71		
Temp. (°C)	4°C			RT		
Conc. (ng/mL)	20	50		20		
N	3			3		
Drug	Butal.	Seco.	Pheno.	Butal.	Seco.	Pheno.
Mean (ng/mL)	18.55	19.03	53.77	20.20	21.93	20.33
%RSD	1.9	10.8	0.9	14.5	2.8	8.4
%Nom	92.8	95.2	107.5	101.0	109.7	101.7

Table 3: Stability Results for Butalbital, Secobarbital and Phenobarbital

Correction Factor

Values reported represent diluted oral fluid. To convert to whole saliva, you must multiply by a factor of 3X.

Conclusions

The ease of use of the Intercept® oral fluid sampling device from OraSure provides an accurate and fast sampling method for many drugs of abuse. The combination of the oral fluid extraction procedure with the analysis speed of the LDTD-MS/MS is an ideal solution in high-throughput drug analysis.

A fast, sensitive and reproducible method for the analysis of barbiturates in oral fluid matrix is achieved using a simple buffer extraction method combined to the speed of analysis of the LDTD-MS/MS with a total sample-to-sample analysis time of **8 seconds**.