

## High-Throughput analysis method for straight chain alkanes Using LDTD-MS/MS

Serge Auger, Alex Birsan, Gregory Blachon and Pierre Picard  
Phytronix Technologies, Québec, Canada,

**Keywords:** High-throughput, alkanes, LDTD

### Introduction

Straight Chain Alkanes are usually analyzed by GC-MS system using Electron Ionization (EI) mode as LC-MS does not properly ionize this type of molecule. To increase the sample throughput, LDTD-MS/MS is used for **10 seconds per sample analysis**.

The Laser Diode Thermal Desorption (LDTD) combined to a mass spectrometer is used to ionize and quantify alkane molecules. The LDTD is a rapid analysis approach in which samples are thermally desorbed. Molecules are channeled, using a carrier gas, to a corona discharge region for ionization prior to detection via a mass spectrometer. Pentacosane and Hexacosane were used as straight chain alkanes to evaluate the ionization process and the development of a quantitative method.

### LDTD-MS/MS System



Figure 1: LDTD system on AB SCIEX 5500 Qtrap Mass Spectrometer

### Method Development

#### Sample preparation

Stock solution of Pentacosane and Hexacosane were dissolved in hexane and the following solutions were prepared:

- 1) Optimization solution of 10 µg/mL in Hexane for both compounds

- 2) Standard curve of Pentacosane using Hexacosane as internal standard.

- 50 ng/mL to 5000 ng/mL
- 4 µL of Sample was added in LazWell plate and evaporated to dryness.

### Hydrocarbon optimization

#### Q1 scan (positive):

Primary mass obtained for Pentacosane (MW: 352.68) was 351 m/z. This mass can be explained by the double bond formation followed by a positive ionization charge: M (-2H) (+H)




Figure 2: Q1 scan mass spectra of Pentacosane in positive mode

#### Product ion scan (pos):

Product ion was then generated with the primary mass, 351 m/z. With a low collision energy (CE:15), typical mass spectra was obtained with loss of 14.




Figure 3: Product ion scan mass spectra of Pentacosane in positive mode

## Dehydrogenation mechanism

A double bond formation during LDTD analysis of Hydrocarbon can be explained by the catalytic mechanism for alkane dehydrogenation reported by Weckhuysen et al<sup>1</sup>. Sample was dried in LazWell plate containing stainless steel sheet. Sample was then vaporized and transferred to the ionization region of the mass spectrometer.

1) Reference : B.M. Weckhuysen, R.A. Shoonheydt, Catalysis Today, 51 (1999) 223-232




Figure 4: Catalytic mechanism reported by Weckhuysen et al.

## LDTD-MS/MS Parameters

### Laser power pattern:

- Increase laser power to 45 % in 3.0 s
- Stay at 45% for 2.0 s
- Decrease laser power to 0 %

Carrier gas flow: 3 L/min (Air)

### MS Parameters

APCI (+)

Scan time: 0.050 s

CE : 22 eV

NC : 3µA

MRM:

- Pentacosane: 351 -> 71
- Hexacosane: 365 -> 71

## Results and Discussion

### Linearity Results

As shown in **Figure 5**, excellent linearity ( $r^2 > 0.99$ ) with no signs of carryover effect is achieved within the quantification range (50 to 5,000ng/mL) for Pentacosane.




Figure 5: Typical standard curve




Figure 6:  
Typical desorption peak:  
A) Blank  
B) STD 50  
C) STD500

## Accuracy and Precision

As shown in **Table 1**, the intra-run accuracy and precision are between 96.8 to 103.8% and 1.3 to 7.8% for Pentacosane.

Drug	S50	S100	S500	S1000	S5000
Nominal conc (ng/mL)	50	100	500	1000	5000
N	3	3	3	3	3
Mean (ng/mL)	49.9	103.8	484.0	990.4	5021.8
%RSD	1.3	4.7	4.7	2.6	7.8
%Nom	99.9	103.8	96.8	99.0	100.4

Table 1: Accuracy and Precision result

## Conclusions

Using the LDTD technology, alkanes can be ionized and quantified according to catalytic mechanism reported by Weckhuysen et al.

A fast analysis can be achieved using LDTD-MS/MS system. This system allows a total sample-to-sample analysis time of **8 seconds**.

For more information about your specific application, visit [www.phytronix.com](http://www.phytronix.com)

Phytronix Technologies  
Parc technologique du Québec métropolitain  
4535, boulevard Wilfrid-Hamel, suite 120, Québec (Qc) Canada G1P 2J7  
[www.phytronix.com](http://www.phytronix.com)