

## LDTD-MS/MS Method validation for Quantitative analysis of small molecule in DBS

Patrice Tremblay, Pierre Picard and Serge Auger

**Keywords:** High-throughput, LDTD, Tandem mass spectrometry, Dextrorphan (DT), FTA-DMPK-C.

### Introduction

A complete method validation was performed for the quantification of Dextrorphan extracted from dry blood spot (DBS). A high-throughput LDTD-MS/MS method was used for quantification of Dextrorphan. The following validation parameter was tested:

- Accuracy and precision of intra and inter-assay
- Matrix selectivity
- Matrix effect
- Recovery
- Stability tests.

### Instrumentation

- Phytronix Technologies LDTD ion source (model T-960);
- TSQ<sup>®</sup> Vantage, Thermo Fisher Scientific.

### LDTD ionization process

The LDTD ion source uses an infrared laser diode to desorb sample that have been dried onto a well of a LazWell™ (96-well plate). The desorbed gas phase molecules are carried into a corona discharge region to undergo APCI, and then they are transferred directly into the mass spectrometer for detection.

### Samples Preparation

#### Blood spiking

- 1 ml Human Blood (Na. Citrate)
- 10µL of stock solution DT in methanol for STD and QC preparation.
- Vortex 0.5 min.

#### Blood spotting

- 15 µl Blood on card.
  - FTA-DMPK-C
- Dry at room temperature for 2h (minimum).

#### Punch extraction

- Add 3mm punch in eppendorf tube (0.5ml)
- 25 µL of NaCl (sat) solution in water
- 50 µL of Internal standard (DT-d3, 5 ng/ml in acetonitrile). Use acetonitrile for Blank.
- Vortex 0.5 min. and centrifuge 2 min. at 14000g.
- Transfer 5.0 µL of upper phase onto LazWell™

### MS Parameters

Mode	APCI (+)
Collision gas	1.5 mTorr (Argon)
S-Lens	30 V
Collision energy	95 V
Scan time	0.020 sec
Needle current	3 µA
Dextrorphan	258.2 → 199.1 amu
Dextrorphan-d3	261.2 → 199.1 amu

### LDTD Parameters

Laser power pattern: 0 to 45% in 3.0 sec.

Carrier gas flow: 3 L/min (Air)

## Results and Discussion

### Calibration Curves

Quantitative determination of Dextrorphan extract for DBS on card FTA-DMPK-C can be achieved over a nominal concentration range of 2 to 500 ng/ml for DT with FTA-DMPK-C card (Figure 1). An excellent linearity is obtained over the concentration range ( $R^2 > 0.99$ ) in three different run.




Fig. 1 Calibration curve of FTA-DMPK-C

	Day 1	Day 2	Day 3
$R^2$	0.9967	0.9985	0.9935
Slope (ratio area / concentration)	0.0157	0.0074	0.0068
y-Abciss	0.0060	0.0025	0.00095

Table 1 Calibration curve parameter.

## Results and Discussion

### Accuracy and Precision (Intra and Inter-assay)

Five levels of QC samples were analyzed in sixuplicate to evaluate the LDTD-MS/MS method accuracy and precision for the intra-assay. The accuracy was evaluated to be within 94.25 and 111.04 % and the precision was within 3.95 and 8.61 % (Table 2)

Intra-assay	LLOQ	QC (Low)	QC (med)	QC (High)	ULOQ
Nom. Conc (ng/ml)	2	8	63	250	500
N	6	6	6	6	6
Mean (ng/ml)	2.2	8.1	59.4	261.9	499.9
RSD (%)	8.6	6.2	4.0	4.0	6.1
%Nom.	111.0	101.0	94.3	104.8	100.0

Table 2 Intra-run accuracy and precision for cards

Three levels of QC samples were analyzed in sixuplicate in three different run to evaluate the LDTD-MS/MS method accuracy and precision for the inter-assay. The accuracy was evaluated to be within 95.16 and 101.76 % and the precision was within 5.38 and 8.99 % (Table 3).

Inter-assay	QC (Low)	QC (med)	QC (High)
Nom. conc. (ng/ml)	8	63	250
N	18	18	18
Mean (ng/ml)	8.1	60.0	254.4
RSD (%)	9.0	6.5	5.4
%Nom.	100.6	95.2	101.8

Table 3 Inter-run accuracy and precision for cards

### Recovery

The percentages of recovery were evaluated for card FTA-DMPK-C.

FTA-DMPK-C	QC (Low)	QC (High)
%Recovery	98.9%	84.0%

Table 3 Recovery of DT at low and high level.

### Matrix selectivity

Six different human blood matrix (2 females and 4 males) was evaluated. The percentage of interference of each blank was evaluated blank peak area compare to the mean peak area value of LLOQ. All blank are lower 20%.

Blank ID	%Interference LLOQ
B1-F	18,8
B2-F	10,0
B3-M	2,0
B4-M	10,9
B5-M	14,1
B6-M	16,7

Table 4 Matrix selectivity evaluation of six different blank

### Matrix effect

Six different matrix were spiked at low QC level and extracted in triplicate. The accuracy was evaluated to be within 91.10 and 104.28 % and the precision was within 6.73 and 14.50 % (Table 5).

Matrix ID / Type	Nom.conc. (ng/ml)	N	Mean (ng/ml)	RSD (%)	%Nom. conc.
M1-F	8	3	8.3	14.5	104.3
M2-F	8	3	7.3	6.7	91.1
M3-M	8	3	7.9	11.2	98.9
M4-M	8	3	8.0	6.8	100.1
M5-M	8	3	8.3	8.3	104.1

Table 5 Matrix effect evaluation

### Stability test result

A pourcentage of deviation from initial value was evaluated for different stability test and a mean value of %RSD was reported (Table 6).

Stability test	Dextrophan	
	QC (Low)	QC (High)
Bench top in blood (94h, 4 °C)		
Mean (%Difference)	7.9	9.6
Precision (Mean %RSD)	6.8	3.3
Dry on FTA-DMPK-C (120h, RT)		
Mean (%Difference)	-0.7	1.3
Precision (Mean %RSD)	9.2	4.3
Extraction solution (88h, 4 °C)		
Mean (%Difference)	11.0	7.6
Precision (Mean %RSD)	7.2	4.9
Dry in LazWell (66h, RT)		
Mean (%Difference)	2.5	--
Precision (Mean %RSD)	13.8	--

Table 6 Stability result

## Conclusions

A complete method validation was performed using dry blood spot (DBS) extract and a LDTD-MS/MS method. LDTD-MS/MS allows high-throughput quantification of Dextrophan from DBS on card FTA-DMPK-C with a sample-to-sample run time of 10 seconds.