

OVERVIEW

Purpose

- High-throughput screening of 16 Benzodiazepines in urine.

Method

- Enzymatic hydrolysis followed by Liquid-Liquid extraction was used for the Benzodiazepines analysis.

Quantification:

- Linearity: $r^2 > 0.995$ over the calibration range (50 to 1000 ng/mL)

- Samples were analyzed with a run time of 6 seconds using LDTD-MS/MS system**

INTRODUCTION

Toxicology laboratories generally use screening methods to obtain a semi-quantitative response for drug samples. Some screening techniques are fast but lack specificity and generate by far too many false positive results. Confirmation of those additional false positive samples is both time and cost consuming. Using mass spectrometry combined with high-throughput LDTD[®] Ion Source enhances specificity at equivalent or better speed. Method assessment is achieved by cross validation with LC-MS/MS, the standard gold method, on the same sample extracts. LC runs were adapted to crude sample preparation by using a 30 minute run in order to reduce ionic suppression. Purified beta- glucuronidase enzymes are used to reduce incubation time to 15 minutes (instead of 1 hour in the original method). Comparison with conventional glucuronidase enzyme incubation is performed in order to validate the obtained results. Complete workflow uses TECAN robotic system with 8 channels liquid handler (figure 1).




Figure 1 Tecan Freedom Evo[®] robot system

LDTD[®] Ionization Source:

The LDTD[®] uses a Laser Diode to produce and control heat on the sample support (Figure 2) which is a 96-well plate. The energy is then transferred through the sample holder to the dry sample which vaporizes prior to being carried by a gas in a corona discharge region. High efficiency protonation and strong resistance to ionic suppression characterize this type of ionization, and is the result of the absence of solvent and mobile phase. This allows for very high throughput capabilities of 6 seconds sample-to-sample analysis time, without carry over.




Figure 2 Schematic of the LDTD ionization source

METHOD

Enzymatic hydrolysis:

- 50 µL patient sample (or standard)
- 15 µL of purified Beta- glucuronidase enzyme
- 25 µL of rapid hydrolysis buffer containing IS solution (200 ng/mL in MeOH)
- Vortex

Incubate 15 minutes at 55°C (131 F)

Liquid Liquid extraction procedure:

- 25 µL of Na₂CO₃ 0.5M pH 10 buffer
- 400 µL Ethyl Acetate/Hexanes 1:1
- Vortex
- Wait for phase separation
- Transfer 4 µL of organic upper layer in LazWell™ plate*
- Dry prior to analysis
- *LazWell™ plate coating:
- 96-well plates for analysis are pre-coated with 5 µL of an EDTA solution (100 µg/mL in MeOH/H₂O/NH₄OH (75/20/5%)) and dried before sample deposition. This coating significantly improves the Temazepam, Lorazepam and Oxazepam signal.

METHOD

LDTD Parameters

- Laser power pattern :
 - Increase laser power to 65 % in 3.0 sec
 - Maintain for 1 sec
 - Decrease laser power to 0 %
- Carrier gas flow (Air) : 3 L/min

MS Parameters

- APCI (+)
- Dwell: 5 msec
- Corona discharge: 3 µA
- DP: 100 V
- MRM mode (see Table 1)

Instrumentation

- LDTD model: S-960
- MS: Sciex 5500 QTrap[®]



Figure 3 LDTD system on Sciex 5500 QTrap[®]

Compound	Q1	Q3	CE (V)
Nordiazepam	271.1	140.2	32
7-Aminoflunitrazepam	284.1	236.0	36
Diazepam	285.0	154.1	32
7-Aminoclonazepam	286.1	222.2	30
Oxazepam	287.0	240.5	30
Estazolam	295.0	205.0	48
Temazepam	301.1	254.6	25
Alprazolam	311.0	274.0	40
Lorazepam	321.0	275.0	23
α-OH-Alprazolam	325.1	204.9	54
2-OH-Ethylflurazepam	333.1	211.2	46
α-OH-Midazolam	342.1	203.0	35
α-OH-Triazolam	359.0	331.0	36
Chlordiazepoxide	300.0	227.0	35
Clonazepam	316.0	214.0	50
Flunitrazepam	314.0	240.0	40
D5-Oxazepam	292.0	245.9	32
D5-Temazepam	306.1	259.6	25
D4-α-OH-Triazolam	363.1	335.0	36

Table 1 MRM method transitions

LC method parameters:

- Column: SilliaChrom[®]SB-C18 (4.6X200 mm)
- Flow rate: 0.5 mL/min
- Mobile Phase A: H₂O/MeOH (90/10)+1% Formic Acid
- Mobile Phase B: H₂O/MeOH (10/90)+1% Formic Acid
- Volume injection: 5 µL
- Ionization mode: Electrospray Ionization (ESI)

Time (min)	Mobile Phase A (%)	Mobile Phase B (%)
0	100	0
10	0	100
15	0	100
16	100	0
30	100	0

Table 2 LC flow gradient

RESULT

Linearity results:

A standard calibration curve (with all 16 drugs) ranging from 50 to 1000 ng/mL has been prepared in blank urine matrix and analyzed in triplicate. All curves have 0.995 coefficients or better. Figure 4 and 5 present typical calibration curves for Oxazepam with LDTD and LC systems respectively.




Figure 4 Oxazepam calibration curve with LDTD




Figure 5 Oxazepam calibration curve with LC

RESULT

LC-MS/MS cross validation:

The most important aspect in a screening method is to provide a Positive flag for all samples that contain targeted drugs. Using both enzyme preparations, no false negative reports were observed using LDTD-MS/MS with the 38 real patient samples tested. During the assay, false positive results were observed on 2 particular samples. A closer look at the integration shows a reduced signal for all internal standards under the acceptance level. Suppression effects are too significant to have an adequate quantitation. Further analysis shows that those 2 samples contain extreme concentrations of opiates (hundreds of µg/mL) causing this effect. The use of deuterated internal standards corrects the quantitation for all the other samples. A threshold level of IS area is used for the identification of overdosed urine sample.




Figure 6 Real samples screening values compared to LC-MS/MS definitive measurement

Purified vs standard β-glucuronidase:

Two β-glucuronidase have been tested. The purified enzyme allows an incubation time of 15 min instead of 60 min with the other one. All real samples have been analysed with both of the enzymes. Figure 7 below compares the amount of several drugs for the same positive patient sample with the two enzymes.




Figure 7 Comparative LDTD-MS/MS Analysis of Sample #100236 for two different types of β-Glucuronidase containing Nordiazepam, Alprazolam, Temazepam, Oxazepam, OH-Alprazolam and 7-Aminoclonazepam

CONCLUSION

- Simultaneous Screening of 16 Benzodiazepines in urine is performed in **6 seconds sample-to-sample** by LDTD-MS/MS
- Good precisions and accuracies are obtained. No carryover was observed
- Cross-validation with LC-MS/MS shows **no false negative** results
- Test of the two β-glucuronidase enzymes gives equivalent results. A fourfold decrease of the analysis time is found using the purified one
- Workflow using the Tecan robotic system allows an analysis capability of 400 samples per hour